

EURONDA SPA  
36030 Montecchio Precalcino (Vicenza) Italy  
via dell'Artigianato, 7  
tel. +39 0445-329811  
fax. +39 0445-865246  
info@euronda.com  
www.euronda.com

DART

MADE IN ITALY



Euronda  
memory  
system



**E-memory®: a small high technology instrument for an easy data saving and management of all the cycles data**

EURONDA memory system is the ideal instrument for the total traceability of the sterilization and test control cycles carried out by the E9 INSPECTION. It attributes automatically a consequential number and archives all the parameters relative to the cycles carried out, even by more E9 INSPECTION autoclaves, being able to recognize their serial numbers. The values are then transferred to a PC (through a USB port) and thanks to its own software, they can then be managed, printed or forwarded by e-mail. Every dental practice can personalize the report with header, graphic logo and personal data. It is also possible to attribute the name of the operator responsible for the sterilization process and to insert the lot number of the tests carried out in order to always obtain the correspondence between each cycle-operator and test-lot number.

**Kit contents:**

- Data logger
- Holder
- UBS cable (3 mt.)
- Instruction manual and warranty (valid for 2 years)
- Software cd rom